

Porter **Sentry**Sedate^{HP}

Portable 4 Cylinder Nitrous Oxide and Oxygen System



Porter MXR-D Digital and Analog Flowmeter



Porter MXR-1 Analog Flowmeter

Porter's newest Nitrous Oxide Sedation System features portability and security with:

- Dual Locking Doors
- Holds 2 Oxygen and 2 Nitrous Oxide "E" Cylinders
- Oxygen/Nitrous Oxide Quick Connect Ports
- Internal and External Storage
- Working Surface
- Durable Locking Wheels

Models

SENTRY-H34AV
MXR-1 Hospital Sentry Sedate with AVS
SENTRY-H44AV
MXR-D Hospital Sentry Sedate with AVS

Dimensions & Weight

Height - 33"
Depth - 24"
Width - 18"
Weight - 96 lbs (add 2 lbs for Porter Digital MXR) - does not includes tanks



www.porterinstrument.com/medical
1-888-723-4001

Sentry Sedate^{HP} Features For Portability, Security And Convenience



The sturdy cabinet design features both internal and external storage.



Sentry Sedate holds two nitrous oxide cylinders and two oxygen cylinders.



A side hook conveniently holds the breathing circuit and hoses.



Large locking wheels allow the unit to be easily moved from room to room.



The dual locks ensure safe, secure cylinder storage.



Oxygen and nitrous oxide quick connect ports bypass cylinders when wall gas is available.